

QuikBow®

STREAMLINING Intraoperative Traction

QuikBow® is a sterile-ready, single-use tension bow for skeletal traction in the operating room.



Indications:

- Pelvic ring fractures
- Acetabular fractures
- Femoral neck and shaft fractures
- Periprosthetic / interprosthetic fractures
- Tibial fractures
- Calcaneal fractures



Spring loaded
QuikCams™
allow for solo
application.

QuikBow® streamlines intraoperative traction while reducing reliance on CSP, controlling the risk of further OR delays, and making the procedures easier with consistent traction.



2520 Columbia St
Vancouver, BC V5Y 3E9


arbutusmedical.com
+1-778-296-2577
info@arbutusmedical.ca

Get QuikBow®



QuikBow®

STREAMLINING Intraoperative Traction



Ask us about our recycling program!



Improve OR efficiency

QuikBow® eliminates costly OR delays caused by underperforming Kirschner bows and replacement equipment wait times



Modern Tension Bow Design

The sterile-ready, single-use tension bow is

- Recyclable
- Never stripped
- Always available




Alleviate CSP Burden

QuikBow® is sterile-ready and single-use, requiring no reprocessing



Reduce Imaging Artifacts

QuikBow® is 99% bio-based polymer, reducing artifacts during imaging



"Intraoperative skeletal traction often makes all the difference in making a hard case easy. Unfortunately, conventional K-bows make what should be a quick, easy step hard and frustrating. I never want to use a conventional K-bow again."

- Allen Nedley MD, OU Health

Get QuikBow®



2520 Columbia St
Vancouver, BC V5Y 3E9

arbutusmedical.com
+1-778-296-2577
info@arbutusmedical.ca

