

Everything you need for **Skeletal Traction** in one procedure pack.



TrakPak® is a complete, single-use skeletal traction system. No specialty beds, no traction frames, nothing to chase down. Just one procedure pack and our surgical power kit.



Eliminate setup scavenger hunt

Reduce pre-procedure setup by 45+ minutes by eliminating the need to chase down traction equipment and supplies.



Standardize the setup workflow

TrakPak® removes tribal knowledge and cross-departmental dependencies from every setup.



Scale without adding Chaos

TrakPak® keeps setup consistent at any volume. No traction beds, no frames, no borrowed equipment to run out of.



**Watch TrakPak®
in Action**

Grab & Go Skeletal Traction

Used in
90+
facilities

What's in the pack?

- 1 SteriGo® power tool drape
- 1 SteriTrak® traction pin driving system 2mm x 310mm
- 1 QuikBow® pin tensioner
- 1 set of Pin Tip Covers
- 4 Blue Towels
- 4 Sticky Drapes
- 2 10.5ml of ChloraPrep w/ Applicator
- 1 Roll of Kerlix
- 1 Traction Rope
- 2 18ga x 1.5 Needle
- 2 22ga x 1.5 Needle
- 2 Xeroform Pads
- 1 20cc Syringe
- 1 #11 Blade
- 1 Skin Marking Pen
- 4 4" x 4" Gauze Pads
- 2 Specialized Adaptors
- 1 Water Bag (20lbs)



The Power Behind TrakPak®

- 1 modified DEWALT drill
- 2 batteries, 1 charger
- 1 pin snips, 1 carry bag



**Watch TrakPak®
in Action**

"Having to leave the room in the middle of the procedure often impedes workflow, elongates the procedure, and risks breaking sterility. The Arbutus Medical kit definitely helps solve this problem."

Greg Pereira MD
Duke University Hospital